

Dimensions: 7 and 2 Series

A. 7 and 2 Series (top view)

B. 2 Series (with doors)

C. 7 Series (with doors)

Height (7 Series)

With Castors	78" (200cm)
Without Castors	74½" (192 cm)

Height (2 Series)


With Castors	30" (76cm)
Without Castors	26½" (67 cm)

Depth (7 and 2 Series)

Wall Mounted <i>without doors</i>	4¾" (13 cm)
Free-Standing <i>without doors</i> <i>including feet that extend in front of the unit</i>	12" (31 cm) 23½" (60 cm)
Wall Mounted <i>with doors</i>	8½" (22 cm)
Free-Standing <i>with doors</i> <i>including feet that extend in front of the unit</i>	15" (39 cm) 23½" (60 cm)

Width (7 Series)

All Units <i>with doors closed</i>	21" (54 cm)
All Units <i>doors fully extended</i>	39½" (100 cm)



Lamp Output:

All readings are measured in milliwatts per square centimeter (mW/cm²).

Approximate median outputs based on historical data.

Model Line	Number of Lamps	UVA MIN / MAX	UVB-BB MIN / MAX	UVB-NB MIN / MAX
2 Series	4		1.0 / 2.0	1.0 / 2.1
	8			1.6 / 2.7
7 Series	4		1.3 / 2.3	1.6 / 2.9
	6			2.3 / 4.3
	8		1.9 / 3.1	2.3 / 4.1
	10	6.97 / 9.45	3.3 / 3.9	3.0 / 5.7

Installation Specifications:

All 7 and 2 Series	Volts	Circuit Rating	Receptacle	Wire Gauge
USA, Canada, Mexico / 60 hertz	110 – 120	15 Amps	household	as locally specified
International Countries / 50 hertz	220 – 240	16 Amps	as locally specified	as locally specified
International Countries / 60 hertz	220 – 240	15 Amps	as locally specified	as locally specified



The Global Leader in Phototherapy

205 West Bement Street P.O. Box 626, Bryan, OH 43506

Phone: 419.636.6304 Fax: 419.636.1739

800.322.8546

e.mail: sales@daavlin.com

www.daavlin.com

© Copyright 2004 by The Daavlin Company. All Rights Reserved.

LIT0040 rev0(12/04)