

# Daavlin

Let's Be Clear  
Daavlin *Phototherapy Solutions*

*7series*



## The Daavlin 7 Series

The 7 Series makes full body phototherapy treatment available in offices with limited space or patient flow.

## The Daavlin 7 Series IS EFFICIENT, ACCURATE AND ELEGANT

The Daavlin 7 Series devices are specifically designed to provide full-body treatment in a small space.

### Efficient

- Our 7 Series X12 utilizes Daavlin's signature double door design to maximize the amount of light that the skin receives and facilitates treatment of all body shapes and sizes.
- Door lamps double body coverage and cut treatment times in half. These doors also conceal and protect the lamps from breakage when the unit is not in use.
- Available lamp types are UVA, Narrow Band UVB or Broad Band UVB.

### Accurate

- Daavlin's optional integrating dosimetry system provides consistent treatment by constantly measuring the energy output in real time, ensuring that you are delivering the correct dose, regardless of factors which may affect treatment times, such as lamp age, room temperature or line voltage.
- Integrating dosimetry automatically compensates for variations in output and delivers the dose you expect, every time.
- Optional DistanceMinder Patient Positioning Sensor ensures optimum placement of your patient with respect to the lamps.
- A security "lock-out" feature prevents unauthorized use of the device. UV blocking eyewear is provided with each unit and protective grid prevents direct contact with the lamps.

### Elegant

Surprisingly slim, all 7 Series devices require just over three square feet of floor space.

Castors and a space-saver wall mounting kit are available as options.



### Daavlin 7 Series Data

|   |             |                 |
|---|-------------|-----------------|
| Treatment Area (no door lamps)                            | 19x70.25 in | (48.3x178.4 cm) |
| Treatment Area (with door lamps)                          | 33x70.25 in | (83.8x178.4 cm) |
| Height  | 74.5 in     | (189.2 cm)      |
| Width (no door lamps)                                     | 20 in       | (50.8 cm)       |
| Width (with door lamps)                                   | 39.5 in     | (100.3 cm)      |
| Depth (including feet)                                    | 24 in       | (61 cm)         |
| Depth (without feet)                                      | 6.5 in      | (16.5 cm)       |
| Weight (no doors)   | 75 lb       | (34 kg)         |
| Weight (with doors)                                       | 106-129 lb  | (48-58.5 kg)    |
| Available in 110-120 V, 60 Hz and 220-240 V, 50 or 60 Hz. |             |                 |

**For further information please contact Daavlin:**

Phone: 1.800.322.8546 Fax: 419.636.1739

International: 001.419.636.6304 Website: [www.daavlin.com](http://www.daavlin.com)

All Daavlin phototherapy units are 510k cleared by the FDA.