

Daavlin

Let's Be Clear
Daavlin Phototherapy Solutions

Mseries



The Daavlin M Series

The Daavlin M Series is a convenient, versatile and technologically advanced system for treatment of photo-responsive diseases of the hands and feet.

The Daavlin M Series MEANS SHORTER TREATMENT TIMES

The Daavlin M Series hand and foot devices are smaller and lighter weight than comparable competitive units. The use of new, high efficiency lamps means that treatment times are dramatically shorter as well.

- The M Series comes with a complement of ten lamps and can be equipped with UVA and / or Narrow Band UVB lamps.
- Combination devices are equipped with six UVA and four NB UVB lamps.
- The lamps used in this device double the power output and effectively eliminate “low energy” spots at the ends of the lamps, as commonly seen in competitive units.

Technology

Daavlin’s optional Flex Dosimetry offers the greatest treatment consistency.

- Doses are entered in Joules for UVA and milliJoules for UVB.
- Conditions that affect output such as lamp age and room temperature are sensed and their effects are automatically integrated into the dose time.

An alternative control option is Daavlin’s Flex Digital Timer.

- This system, in which doses are entered in minutes and seconds may be attractive if a meter for measuring the output of the device is already available.



Safety Features

- A slide-on hood lined with spectral aluminum enhances efficiency by reflecting ambient UV energy back into the treatment area.
- Two lamps in the hood enable simultaneous treatment of both aspects of the hands and or feet.
- An optional table accommodates two units and provides the flexibility of simultaneous treatment of the hands and feet.
- Equipped with castors, the table makes the system simple to store and frees up valuable office space when the device is not in use.
- The hood is removable which allows the device to be used as a panel unit when turned on its side.

Technology

UV transmitting acrylic provides a comfortable resting place for the hands or feet while preventing fingers or toes from direct contact with the lamps, as seen with competitive devices.

Daavlin M Series Data

Treatment Area	16.5x17.75 in	(41.9x45.1 cm)
Height (with hood)	13 in	(33 cm)
Width	20.25 in	(51.4 cm)
Depth	22 in	(55.9 cm)
Weight	48 lbs	(21.8 kg)
Table Dimensions:	26"W x 26"D x 33.5"H (66 cm x 66 cm x 85 cm)	
Available in 110-120 V, 60 Hz and 220-240 V, 50/60 Hz.		

For further information please contact Daavlin:

Phone: 1.800.322.8546 Fax: 419.636.1739

International: 001.419.636.6304 Website: www.daavlin.com

All Daavlin phototherapy units are 510k cleared by the FDA.